




QUOTIENT

ALBAcyte® A₂ Cells REAGENT RED CELLS



INTERPRETATION OF LABEL SYMBOLS

	Batch code
	Use by (YYYY-MM-DD)
	Storage temperature limitation (2-8 °C)
	<i>In vitro</i> diagnostic medical device
	Product code
	Consult instructions for use
	Manufacturer

INTENDED PURPOSE

These reagent red cells are for the ABO reverse grouping of patient or donor serum/plasma and may be used for the control of serological tests.

INTRODUCTION

ABO blood grouping is generally performed by testing red cells with anti-A and anti-B (many laboratories also test with anti-A,B). A₁, A₂, B and O cells should be used to control each batch of tests. Confirmation of the red cell group can be provided by simultaneously performing a reverse or serum group i.e. testing the donor or recipients serum/plasma with reagent red blood cells of groups A₁ and B to detect anti-A and anti-B. Group A₂ reagent red blood cells can be used to identify anti-A₁ in the serum of group A people. Group O cells can be used to identify agglutination due to non-ABO agglutinins.

REAGENT DESCRIPTION

These reagent red cells are presented as a 2-3% suspension of washed red cells in Modified Alsever's Solution. The cells are group A₂. The preservative solution has been specially formulated to preserve red cell integrity and antigenicity and contains the following components - trisodium citrate, citric acid, dextrose, inosine, neomycin sulphate (0.103g/l) and chloramphenicol (0.349g/l).

The volume delivered by the reagent dropper bottle is approximately 40µl; bearing this in mind, care should be taken to ensure that appropriate serum: cell ratios are maintained in all test systems.

This reagent complies with the requirements of Directive 98/79/EC on *in vitro* Diagnostic Medical Devices and the recommendations contained in the Guidelines for Blood Transfusion Services in the United Kingdom.

PRECAUTIONS FOR USE AND DISPOSAL

Source material from which this product is derived was found non reactive for HBsAg, Anti-HIV 1/2 and Anti-HCV.

No known test method can offer assurances that products derived from human blood will not transmit infectious disease; therefore appropriate care should be taken in the use and disposal of this product.

Chloramphenicol is classified as a carcinogen and neomycin sulphate is classified as an irritant.

This reagent is for *in vitro* professional use only.

These reagent red cells may be used directly from the vial or may be washed and resuspended before use to 2-3% in PBS or 1.5-2% in LISS. Reagent red cells treated in this way must be discarded within 24 hours of preparation. Transfer of these reagent red cells for long term storage in another container is not recommended.

Furthermore, when the user changes the reagent in any way, e.g. the preparation of LISS cell suspensions, the user is responsible for assuring the strength of red cell suspension, the quality of PBS or LISS used and the generation and storage of relevant documentation

STORAGE CONDITIONS

The reagent should be stored at 2-8 °C. Do not freeze. Do not use if obviously discoloured or haemolysed. Do not use beyond the notified expiry date.

SPECIMEN COLLECTION AND PREPARATION

Specimens should be collected by aseptic technique. The specimen should be tested as soon as possible after collection. If testing is delayed, the specimen should be stored at 2-8 °C for a maximum of 48 hours. Blood specimens exhibiting gross haemolysis or contamination should not be used.

TEST PROCEDURES

No specific test procedures are recommended. Users are advised to carefully validate procedures and confirm reagent suitability before use.

PERFORMANCE LIMITATIONS

The presence of irregular antibodies in the serum/plasma of a patient/donor may cause unexpected agglutination of these reagent red cells.

Some loss of antigenic expression may occur during the stated shelf life. Since this loss is partly determined by characteristics of individual blood donations or donors which cannot be predicted or controlled, the recommended conditions of storage and use must be rigidly applied.

False positive or false negative results can occur due to contamination of test materials, improper reaction temperature, improper storage of materials, omission of test reagents and certain disease states.

A₂ cells have been characterised using the Anti-A₁ lectin *Dolichos biflorus*. This lectin in its undiluted form has Anti-A specificity and requires appropriate dilution to react directly with A₁ and A₁B red blood cells yet fail to react with A₂ and A₂B red blood cells. A small number of red blood cells may give strong reactions with

Anti-A and weak reactions with Anti-A₁. The status of these red blood cells may be regarded as A₂ a subgroup with phenotype characteristics in between A₁ and A₂. It should be noted that some A₂ and A₂B cells may react weakly with Anti-A₁ if incubated too long. Anti-A₁ lectin would be expected to react 2+4+ with A₁ and A₁B red blood cells.

SPECIFIC PERFORMANCE CHARACTERISTICS

The reagent red cells have been shown to have a negative direct antiglobulin test, indicating that no human IgG or C3 complement components are detectable on the cell surface.

In performance evaluation studies (data on file at Alba Bioscience Limited), Z406 was tested against random plasma samples and ABO antisera. The performance of Z406 is summarised as positive and negative percentage agreement.

Positive percentage agreement was 100% and negative percentage agreement was 100%.

The test outcomes were 100% in concordance with the expected test outcome based on the antibody contained in the plasma or reagent.

DATE OF ISSUE

2022-07

For further information or advice please contact your local distributor.



Emergo Europe B.V.
Prinsessegracht 20
2514 AP The Hague
The Netherlands



Alba Bioscience Limited
James Hamilton Way,
Penicuik, EH26 0BF
UK

Tel No: +44 (0) 131 357 3333
Fax No: +44 (0) 131 445 7125
E-Mail: customer.serviceEU@quotientbd.com

© Alba Bioscience 2022

Z406PI/08

INSTRUCTIONS FOR USE

Instructions for Use are available on Quotient website at www.quotientbd.com or can be requested from your Local Distributor, by providing the relevant Product Code stated on product labels and Instructions for Use supplied with the product.

NOTICE D'UTILISATION

La notice d'utilisation est disponible sur le site Internet de Quotient à l'adresse www.quotientbd.com ou peut être demandée auprès de votre distributeur local, en précisant le code produit correspondant indiqué sur les étiquettes du produit et la notice d'utilisation fournie avec le produit.

ISTRUZIONI PER L'USO

Le istruzioni per l'uso sono disponibili sul sito Web di Quotient, www.quotientbd.com, o possono essere richieste al proprio distributore locale fornendo il codice del prodotto indicato sulle etichette del prodotto e nelle istruzioni per l'uso fornite con il prodotto.

GEBRAUCHSANWEISUNG

Die Gebrauchsanweisung ist auch auf der Quotient Internetseite unter www.eu.quotientbd.com/ifu erhältlich oder kann bei Ihrem zuständigen Vertriebspartner angefordert werden. Hierfür geben Sie bitte die jeweilige Artikel-Nummer an, die sich auf den Etiketten und Gebrauchsanweisungen befindet, die mit dem Produkt geliefert werden.