



1 x 10 mL or 3 x 10 mL

ALBAcyte® IgG Sensitized Red Blood Cells (Z441U) ALBAcyte® IgG Sensitized Red Blood Cells II (Z443U)

Control the indirect (IAT) and direct (DAT) antiglobulin tests. Available in two formulations to allow you to select the strength of reactivity that aligns with your laboratory's preference.



4 x 10 mL

ALBAsure™ QC Kit (Z481U)

Simplify daily quality control of manual ABO, RhD, and antibody detection reagents.

QC Kit Cell 1 AB rr red cells (2-3%)

QC Kit Cell 2 O R1r red cells (2-3%)

QC Kit Antibody (x2) Antibody blend (Anti-A, -B, -D, -c)



1 x 3 mL

ALBAcyte® C3 Coated Red Blood Cells (Z482U)

Confirm the presence or absence of Anti-C3 reactivity of the Anti-Human Globulin Reagent and use as indicator cells in AHG testing.



1 x 10 mL or 10 x 10 mL

ALBAcheck® Monoclonal Control (Z271U)

Use as a negative control for select Alba by Quotient monoclonal blood grouping reagents.



TO REQUEST A QUOTE OR SAMPLES, SCAN THIS CODE.

Visit us at quotientbd.com or email us at Customer.ServiceUS@quotientbd.com

