

QUALITY BLOOD BANK REAGENTS

The Alba by Quotient manual reagent product line delivers solid value while maintaining quality. Join Quotient's satisfied customers and see for yourself the difference that Alba can make to your laboratory.

Building on over 30 years in transfusion diagnostics, Alba by Quotient supplies reagents to over 5,000 labs around the globe.

ABO/RHD ANTISERA

PRODUCT NAME	PACK SIZE	PRODUCT CODE
ALBAclone® Anti-A	1 x 10 ml or 10 x 10 ml	Z001U-1 or Z001U-10
ALBAclone® Anti-B	1 x 10 ml or 10 x 10 ml	Z015U-1 or Z015U-10
ALBAclone® Anti-A,B	1 x 10 ml or 10 x 10 ml	Z023U-1 or Z023U-10
ALBAlect™ Anti-A ₁ lectin	1 x 5 ml	Z241U
ALBAclone® Anti-D fusion (IgM/IgG)	1 x 10 ml or 10 x 10 ml	Z043U-1 or Z043U-10
ALBAclone® Anti-D blend (IgM/IgG)	1 x 10 ml or 10 x 10 ml	Z041U-1 or Z041U-10
ALBAclone® Anti-D alpha (IgM)	10 x 10 ml	Z031U
ALBAcheck® BGS Monoclonal Control	1 x 10 ml or 10 x 10 ml	Z271U-1 or Z271U-10

REAGENT RED BLOOD CELLS

> REVERSE GROUPING

PRODUCT NAME	PACK SIZE	PRODUCT CODE
ALBAcyte® A ₁ Cells	1 x 10 ml	Z401U
ALBAcyte® A ₂ Cells	1 x 10 ml	Z406U
ALBAcyte® B Cells	1 x 10 ml	Z411U
ALBAcyte® Orr Cells	1 x 10 ml	Z421U

> ANTIBODY DETECTION & IDENTIFICATION

PRODUCT NAME	PACK SIZE	PRODUCT CODE
ALBAcyte® Antibody Screen (3-cell)	3 x 10 ml	Z451U
ALBAcyte® Antibody Screen (2-cell)	2 x 10 ml	Z454U
ALBAcyte® Expanded Rh Negative Antibody Screen	4 x 5 ml	Z464U
ALBAcyte® Antibody ID Panel (10-cell)	10 x 5 ml	Z471U
ALBAcyte® Antibody ID Panel (10-cell) papain-treated	10 x 5 ml	Z472U
ALBAcyte® Antibody ID Panel (16-cell)	16 x 5 ml	Z473U

> QUALITY CONTROL

PRODUCT NAME	PACK SIZE	PRODUCT CODE
ALBAcyte® IgG Sensitized Cells (weak)	1 x 10 ml or 3 x 10 ml	Z441U-1 or Z441U-3
ALBAcyte® IgG Sensitized Cells II	1 x 10 ml or 3 x 10 ml	Z443U-1 or Z443U-3
ALBAsure® QC Kit	4 x 10 ml	Z481U
ALBAcyte® C3 Coated Red Blood Cells	1 x 3 ml	Z482U

ANTI-HUMAN GLOBULIN

PRODUCT NAME	PACK SIZE	PRODUCT CODE
Anti-IgG,-C3d; Polyspecific (Rabbit/Murine Monoclonal) Clear	1 x 10 ml or 10 x 10 ml	Z351U-1 or Z351U-10
Anti-IgG,-C3d; Polyspecific (Rabbit/Murine Monoclonal) Green	1 x 10 ml or 10 x 10 ml	Z350U-1 or Z350U-10
Anti-IgG (Rabbit) Clear	1 x 10 ml or 10 x 10 ml	Z356U-1 or Z356U-10
Anti-IgG (Rabbit) Green	1 x 10 ml or 10 x 10 ml	Z357U-1 or Z357U-10
ALBAclone® Anti-C3d (Murine Monoclonal) (Detects C3b and/or C3d)	1 x 3 ml	Z360U

ANTISERA FOR EXTENDED PHENOTYPING

PRODUCT NAME	TECHNIQUE*	PACK SIZE	PRODUCT CODE
ALBAclone® Anti-C	5-15 minutes	1 x 5 ml or 3 x 5 ml	Z064U-1 or Z064U-3
ALBAclone® Anti-ċ	5-15 minutes	1 x 5 ml or 3 x 5 ml	Z083U-1 or Z083U-3
ALBAclone® Anti-E	5-15 minutes	1 x 5 ml or 3 x 5 ml	Z073U-1 or Z073U-3
ALBAclone® Anti-ē	5-15 minutes	1 x 5 ml or 3 x 5 ml	Z096U-1 or Z096U-3
ALBAclone® Anti-C ^w	5-15 minutes	1 x 2 ml	Z106U
ALBAclone® Anti-K	5-15 minutes	1 x 5 ml	Z132U
ALBAclone® Anti-Fy ^a	IAT, 15-30 minutes	1 x 5 ml	Z152U
ALBAclone® Anti-Fy ^b	Immediate spin	1 x 5 ml	Z154U
ALBAclone® Anti-Jk ^a	Immediate spin	1 x 2 ml	Z162U
ALBAclone® Anti-Jk ^a	Immediate spin	1 x 5 ml or 3 x 5 ml	Z164U-1 or Z164U-3
ALBAclone® Anti-Jk ^b	Immediate spin	1 x 2 ml	Z166U
ALBAclone® Anti-Jk ^b	Immediate spin	1 x 5 ml or 3 x 5 ml	Z168U-1 or Z168U-3
ALBAclone® Anti-M	5 minutes	1 x 5 ml	Z171U
ALBAclone® Anti-N	5 minutes	1 x 5 ml	Z176U
ALBAclone® Anti-S	IAT, 15-30 minutes	1 x 5 ml	Z182U
ALBAclone® Anti-š	IAT, 15-30 minutes	1 x 5 ml	Z187U
ALBAclone® Anti-P1	Immediate spin	1 x 5 ml	Z202U
ALBAclone® Anti-Le ^a	15 minutes	1 x 5 ml	Z212U
ALBAclone® Anti-Le ^b	15 minutes	1 x 5 ml	Z217U
ALBAclone® Anti-Lu ^b	15 minutes	1 x 2 ml	Z223U
ALBAsera® Anti-Wr ^a (polyclonal)	IAT, 15-45 minutes	1 x 2 ml	Z231U

*For reference only — refer to IFU for appropriate technique.

ENHANCEMENT MEDIA

PRODUCT NAME	PACK SIZE	PRODUCT CODE
ALBAhance™ Bovine Serum Albumin (BSA) 22%	10 x 10 ml	Z305U
ALBAhance™ PEG	10 x 10 ml	Z312U
ALBAhance™ LISS Additive	10 x 10 ml	Z333U
ALBAzyme™ Papain Solution Kit (Papain Solution & Enzyme Control)	1 x 3 ml & 1 x 5 ml	Z317U

SPECIALTY PRODUCTS

PRODUCT NAME	PACK SIZE	PRODUCT CODE
ALBAcheck® Competency Testing Kit	10 x 2 ml	Z276
Quotient Reagent Vial Rack (15-vial)	----	Z891U
Quotient Reagent Vial Rack (18-vial)	----	Z892U

All products meet relevant FDA potency standards.



WHY CHOOSE QUOTIENT?



Expert technical support from experienced blood bankers



Samples and protocols for reagent evaluations available



Dedicated customer portal



Customer service you can rely on



Customizable cost savings analysis available

Alba by Quotient's customers benefit from high-quality products that support and complement operational processes.



TO REQUEST A QUOTE OR SAMPLES, SCAN THIS CODE.

Visit us at quotientbd.com or email us at Customer.ServiceUS@quotientbd.com



About Quotient

Building on over 30 years of experience in transfusion diagnostics, Quotient is a commercial-stage diagnostics company committed to delivering solutions that reshape the way diagnostics is practiced. Quotient has two key product portfolios, Alba by Quotient™ and MosaiQ™ by Quotient. The Alba by Quotient product range is comprised of high-quality reagents used to keep laboratories around the world running efficiently. They are sold directly to laboratories or through our OEM partners.

MosaiQ, Quotient's proprietary multiplex microarray technology, offers a fully automated, consolidated testing platform, allowing for multiple tests across different modalities. MosaiQ is designed to be a game-changing solution, which Quotient believes will increase efficiencies, improve clinical practice, deliver significant workflow improvements, and operational cost savings to laboratories around the world.

©Alba Bioscience Ltd, 2022. The Quotient logo, Alba by Quotient, MosaiQ, ALBAclone, ALBAcyte, ALBAlect, ALBAsure, ALBAhance, ALBAzyme, and ALBAcheck are trademarks or registered trademarks of Quotient Ltd or its subsidiaries in various jurisdictions. Menus and capabilities are subject to change. Not all methods may be available in all territories. Subject to regulatory clearance.