



Patient/Sample Details				
Name:	DOB:	Sample Identification:	ABO Rh:	
Diagnosis/History:	Hospital:	DAT:	Anti-IgG:	Anti-C3

Lot No:	V242210	
Expiry Date:	2022.01.17	U.S. License 1807

Cell #	Rh-hr	Donor	Rh-hr								Kell					Duffy		Kidd		Lewis		MNS				P	Lutheran		Additional Antigens			TEST RESULTS				Cell #		
			D	C	E	c	e	f	*V	C ^w	K	k	Kp ^a	Kp ^b	*Js ^a	*Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	M	N	S	s	P1	Lu ^a	Lu ^b	Xg ^a	Wr ^a	Special Types						
1	R ₁ ^w R ₁	2768020030257	+	+	0	0	+	0	NT	+	0	+	0	+	0	+	0	+	+	0	0	+	+	0	+	0	0	0	+	*+	0							1
2	R ₂ R ₂	6803031119013	+	0	+	+	0	0	NT	0	0	+	0	+	0	+	+	+	0	+	+	0	0	+	0	+	+	0	+	*+	*0							2
3	rr	3256732	0	0	0	+	+	+	0	0	+	+	0	+	0	+	+	0	+	0	0	+	+	0	0	+	+	0	+	+	0							3
Patient Cells																																						

Notes:

- All cells are DAT negative.
- NT = Not tested.
- The f antigen status has been determined presumptively based on Rh-hr phenotype

* Indicates those antigens whose presence or absence may have been determined using only a single example of a specific antibody.

Quotient, 301 South State Street, S-204, Newtown, PA 18940, USA
 Product Technical Support Tel: 1-888-228-1990 Web: www.quotientbd.com

Conclusions/Further Testing Required:	
Signature/Initials:	Date: