

ALBAcyte® REAGENT RED BLOOD CELLS FOR DETECTION OF UNEXPECTED ANTIBODIES

Antibody Screen (2-Cell) - REF Z454U



Patient/Sample Details				
Name:	DOB:	Sample Identification:	ABO Rh:	
Diagnosis/History:	Hospital:	DAT:	Anti-IgG:	Anti-C3

Lot No:	V250042	
Expiry Date:	2022.08.29	U.S. License 1807

Cell #	Rh-hr	Donor	Rh-hr								Kell					Duffy		Kidd		Lewis		MNS				P	Lutheran		Additional Antigens			TEST RESULTS				Cell #			
			D	C	E	c	e	f	*V	C ^w	K	k	Kp ^a	Kp ^b	*Js ^a	*Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	M	N	S	s	P1	Lu ^a	Lu ^b	Xg ^a	Wr ^a	Special Types							
1	R ₁ R ₁	4231321	+	+	0	0	+	0	0	0	+	+	0	+	0	+	+	0	0	+	0	+	+	0	+	0	+	0	+	0	0								1
2	R ₂ R ₂	11215993	+	0	+	+	0	0	NT	0	0	+	0	+	0	+	0	+	+	0	+	0	+	+	0	+	+	0	+	+	0								2
Patient Cells																																							

Notes:

1. All cells are DAT negative.
 2. NT = Not tested.
 3. The f antigen status has been determined presumptively.
- * Indicates those antigens whose presence or absence may have been determined using only a single example of a specific antibody.

<p>Conclusions/Further Testing Required:</p>
<p>Signature/Initials: _____ Date: _____</p>

Quotient, 301 South State Street, S-204, Newtown, PA 18940, USA
 Product Technical Support Tel: 1-888-228-1990 Web: www.quotientbd.com

Z454UPI/02