



Patient/Sample Details				
Name:	DOB:	Sample Identification:	ABO Rh:	
Diagnosis/History:	Hospital:	DAT:	Anti-IgG:	Anti-C3

Lot No:	V254361	U.S. License 1807
Expiry Date:	2022.12.19	

Cell #	Rh-hr	Donor	Rh-hr							Kell				Duffy		Kidd		Lewis		MNS				P	Lutheran		Additional Antigens			TEST RESULTS				Cell #					
			D	C	E	c	e	f	V	C ^w	K	k	Kp ^a	Kp ^b	*Js ^a	*Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	M	N	S	s	P1	Lu ^a	Lu ^b	Xg ^a	Wr ^a	Special Types							
1	R ₁ ^w R ₁	2768020946517	+	+	0	0	+	0	NT	+	0	+	0	+	0	+	+	+	+	+	+	0	+	+	+	+	+	+	0	*+	*0								1
2	R ₂ R ₂	2768020042889	+	0	+	+	0	0	NT	0	0	+	0	+	0	+	+	0	+	0	0	0	0	+	0	+	0	0	+	*+	*0								2
3	rr	6901082019015	0	0	0	+	+	+	NT	0	+	+	0	+	0	+	0	+	0	+	+	+	+	0	+	0	+	0	+	0	0								3
Patient Cells																																							

Notes:

- All cells are DAT negative.
- NT = Not tested.
- The f antigen status has been determined presumptively based on Rh-hr phenotype

* Indicates those antigens whose presence or absence may have been determined using only a single example of a specific antibody.

Product Technical Support Tel: 1-888-228-1990 Web: www.quotientbd.com

Conclusions/Further Testing Required:
Signature/Initials: _____ Date: _____