

**ALBAcyte® REAGENT RED BLOOD CELLS FOR USE IN IDENTIFICATION OF UNEXPECTED ANTIBODIES**

**Expanded Rh Negative Antibody Screen - REF Z464U**



Patient/Sample Details				
Name:	DOB:	Sample Identification:	ABO Rh:	
Diagnosis/History:	Hospital:	DAT:	Anti-IgG:	Anti-C3

<b>Lot No:</b>	<b>V258670</b>	<b>U.S. License 1807</b>
<b>Expiry Date:</b>	<b>2023.04.24</b>	

Cell #	Rh-hr	Donor	Rh-hr								Kell				Duffy		Kidd		Lewis		MNS				P	Lutheran		Additional Antigens			TEST RESULTS				Cell #				
			D	C	E	c	e	f	*V	C <sup>w</sup>	K	k	Kp <sup>a</sup>	Kp <sup>b</sup>	*Js <sup>a</sup>	*Js <sup>b</sup>	Fy <sup>a</sup>	Fy <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	Le <sup>a</sup>	Le <sup>b</sup>	M	N	S	s	P1	Lu <sup>a</sup>	Lu <sup>b</sup>	Xg <sup>a</sup>	Wr <sup>a</sup>	Special Types							
1	r'r	1231207	0	+	0	+	+	+	NT	0	+	+	0	+	0	+	+	0	0	+	+	0	+	0	+	+	0	0	+	0	0								1
2	r'r	1098873	0	0	+	+	+	+	0	0	0	+	0	+	0	+	0	+	+	0	0	+	+	0	+	0	+	0	+	+	0								2
3	rr	1179-068-0	0	0	0	+	+	+	NT	0	0	+	0	+	0	+	0	0	+	+	0	+	+	+	0	+	+	0	+	+	0								3
4	R <sub>0</sub> R <sub>0</sub>	2389780	+	0	0	+	+	+	NT	0	0	+	0	+	0	+	0	0	+	+	0	+	0	+	+	+	+	0	+	0	0								4
Patient Cells																																							

**Notes:**

- All cells are DAT negative.
- NT = Not tested.
- The f antigen status has been determined presumptively.
- Cell #4 is designated as R<sub>0</sub>R<sub>0</sub> based on the statistical probability for the source donor population but a genotype of R<sub>0</sub>r cannot be excluded.

\* Indicates those antigens whose presence or absence may have been determined using only a single example of a specific antibody.

<b>Conclusions/Further Testing Required:</b>	
<b>Signature/Initials:</b>	<b>Date:</b>