

ALBAcyte® REAGENT RED BLOOD CELLS FOR USE IN IDENTIFICATION OF UNEXPECTED ANTIBODIES

REF Z471U & Z472U



| Patient/Sample Details | | | | |
|------------------------|-----------|------------------------|-----------|----------|
| Name: | D.O.B: | Sample Identification: | ABO Rh: | |
| Diagnosis History: | Hospital: | DAT: | Anti-IgG: | Anti-C3: |

| | | | |
|--------------|-------------------|-------------|------|
| Lot No: | V260421 & V260423 | | |
| Expiry Date: | 2023.06.19 | U.S License | 1807 |

| Cell # | Rh-hr | Donor | Rh-hr | | | | | | | Kell | | | | Duffy | | Kidd | | Lewis | | MNS | | | | P | Lutheran | | Additional Antigens | | | TEST RESULTS | | | | Cell # | | | | | | |
|---------------|--|------------|-------|---|---|---|---|---|----|----------------|---|---|-----------------|-----------------|------------------|------------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|---|---|----------|---|---------------------|-----------------|-----------------|-----------------|-----------------|---------------|--|--------|--|--|--|--|---|----|
| | | | D | C | E | c | e | f | *V | C ^w | K | k | Kp ^a | Kp ^b | *Js ^a | *Js ^b | Fy ^a | Fy ^b | Jk ^a | Jk ^b | Le ^a | Le ^b | M | N | S | s | P1 | Lu ^a | Lu ^b | Xg ^a | Wt ^a | Special Types | | | | | | | | |
| 1 | R ₁ ^w R ₁ | 1685553 | + | + | 0 | 0 | + | 0 | NT | + | 0 | + | 0 | + | 0 | + | + | 0 | 0 | + | + | + | 0 | + | 0 | + | + | + | + | + | 0 | 0 | | | | | | | 1 | |
| 2 | R ₁ R ₁ | 6232863 | + | + | 0 | 0 | + | 0 | NT | 0 | + | + | 0 | + | 0 | + | + | + | 0 | 0 | + | + | + | + | + | 0 | + | + | *0 | 0 | | | | | | | | | 2 | |
| 3 | R ₂ R ₂ | 1352908 | + | 0 | + | + | 0 | 0 | 0 | 0 | 0 | + | 0 | + | 0 | + | 0 | + | 0 | + | 0 | + | + | + | + | + | 0 | + | + | 0 | 0 | | | | | | | | | 3 |
| 4 | r'r | 1231207 | 0 | + | 0 | + | + | + | NT | 0 | + | + | 0 | + | 0 | + | + | 0 | + | + | 0 | + | + | + | 0 | 0 | + | + | 0 | 0 | | | | | | | | | | 4 |
| 5 | r''r | 533687 | 0 | 0 | + | + | + | + | NT | 0 | 0 | + | 0 | + | 0 | + | + | + | 0 | + | + | + | 0 | + | 0 | 0 | + | + | *+ | 0 | | | | | | | | | | 5 |
| 6 | R ₀ r | 1007-333-6 | + | 0 | 0 | + | + | + | NT | 0 | 0 | + | 0 | + | 0 | + | + | + | 0 | + | 0 | + | + | + | 0 | + | 0 | + | *0 | 0 | | | | | | | | | | 6 |
| 7 | rr | 4211534 | 0 | 0 | 0 | + | + | + | 0 | 0 | + | + | 0 | + | 0 | + | + | 0 | + | 0 | + | + | + | 0 | + | + | + | + | + | 0 | | | | | | | | | | 7 |
| 8 | rr | 11709656 | 0 | 0 | 0 | + | + | + | NT | 0 | 0 | + | 0 | + | 0 | + | + | 0 | + | + | 0 | + | + | + | +w | 0 | + | *+ | 0 | | | | | | | | | | | 8 |
| 9 | rr | 10492320 | 0 | 0 | 0 | + | + | + | 0 | 0 | 0 | + | 0 | + | 0 | 0 | + | + | 0 | + | + | 0 | 0 | + | + | 0 | + | + | + | 0 | | | | | | | | | | 9 |
| 10 | rr | 1580417 | 0 | 0 | 0 | + | + | + | NT | 0 | 0 | + | + | 0 | + | + | + | 0 | + | 0 | + | + | + | + | + | 0 | + | + | + | 0 | | | | | | | | | | 10 |
| Patient Cells | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

- Notes:**
- All cells are DAT negative.
 - M, N, S, s, Fy^a and Fy^b antigens are destroyed or depressed following enzyme treatment.
 - 's' = strong reactor, 'w' = weak reactor, NT = not tested
 - The f antigen status may have been determined presumptively based on the Rh-hr phenotype.
- * Indicates those antigens whose presence or absence may have been determined using only a single example of a specific antibody.

| |
|---------------------------------------|
| Conclusions/Further Testing Required: |
| Signature/Initials: _____ Date: _____ |