

Expected Test Results	Reagents					Reagent Red Cells									QC Kit lot Number:		
	Anti-A	Anti-B	Anti-AB	Anti-D	Rh Control	A1	A2	B	O	Screen Cell 1		Screen Cell 2		Screen Cell 3		Expiration Date:	
										37°	IAT	37°	IAT	37°	IAT		
QC Cell 1 – AB rr	1-4+	1-4+	1-4+	0	0												
QC Cell 2 – O R1r	0	0	0	1-4+	0												
QC Antibody						1-4+	1-4+	1-4+	0	0-2+	1-4+	0-2+	1-4+	0-2+	1-4+		

  

Daily Test	Reagents					Reagent Red Cells									Additional Reagents			
	Anti-A	Anti-B	Anti-AB	Anti-D	Rh Control	A1	A2	B	O	Screen Cell 1		Screen Cell 2		Screen Cell 3		AHG	Enhancement Media	IgG Sensitized Cells
										37°	IAT	37°	IAT	37°	IAT			
Lot No:																		
Expiry Date:																		
QC Cell 1 – AB rr																		
QC Cell 2 – O R1r																		
QC Antibody																		
Satisfactory ✓																		
Unsatisfactory*																		

  

Lot No:																		
Expiry Date:																		
QC Cell 1 – AB rr																		
QC Cell 2 – O R1r																		
QC Antibody																		
Satisfactory ✓																		
Unsatisfactory*																		

  

Lot No:																		
Expiry Date:																		
QC Cell 1 – AB rr																		
QC Cell 2 – O R1r																		
QC Antibody																		
Satisfactory ✓																		
Unsatisfactory*																		

  

Lot No:																		
Expiry Date:																		
QC Cell 1 – AB rr																		
QC Cell 2 – O R1r																		
QC Antibody																		
Satisfactory ✓																		
Unsatisfactory*																		